



03.09

## Small animals

### New accessories for the Video-Compact-System



#### Silicon sheath for the Video-Otoscop.

Soft and smooth for sensitive animals semi disposable and easy clean for professional ear care.

### Special ear speculum

- *the combination between latest video otoscopy and conventional otoscopy*



Ear speculum with open portal for foreign body removal. With silicon hose to protect the ear canal.



## Small animals



Flexible Bronchoscope, silicon grip  
– ergonomical handling –

### **FLB-39-80**

**Broncho Fiberscope,**  
OD:3,9mm, WL:70cm  
with ergonomic silicon grip,  
instrument channel:1,4 mm,  
180°up/100°down angle of  
vision: 95°, standard light cable  
connectors, complete with  
suction connector, density  
tester, cleaning brush and  
transportation case,

compatible

### **B12-810**

**Flexible Biopsy Forceps**

---

## Small animals

Special Micro-Endoscope with  
working and irrigation canal  
for urethra and kidney examination.



### **MI-23-1000**

**Mirco Fiberscope,**  
OD:2,3 mm, WL: 100 cm  
working channel: 1,2/0, 5mm,  
90° deflecting, field of vision:70°  
image bundle: 6000

compatible

### **B10-710**

**Flexible Biopsy Forceps**

## Small animals

### Micro rhinoscope semiflex 2,3 mm, WL 20 cm



high definition fiber optic, 0°,  
rigid and semiflexible,  
working channel: 1,1 mm,  
DIN ocular flexible with standard,  
light cable connectors

---


## Small animals

### Mirror system with integr. cold light fibers

for post rhinoscopy in cats and dogs



**cold light system**



Until now examination of retro pharynx was only possible for clinics an hospitals with good endoscopic equipment. The newly developped deflecting mirror system allows the average practice to perform retro pharyngoscopy. Grasping of foreign bodies or taking biopsies under visual controle are possible. The universal connector of the cold light carrier allows the use of standard light sources and cables or even batterie operated light handles.



## Recorder notebook

(highend) set, "high resolution"

digital image capture & video recording, incl.  
software (USB) and completely installed  
database, "ready to go", compatible to win/XP,  
including video adaptors, compatible to  
network systems.

\*\*\* compatible to conventional endoscopic  
video cameras

---